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92-202 Łódź, Al. Marszałka J. Piłsudskiego 84
tel. (+48 42) 636 92 11, (+48 42) 636 55 72, (+48 42) 674 64 14 wew. 320, fax (+48 42) 674 81 24
zj@ibprs.pl
NIP: 525-000-82-64 REGON: 000053835-00026

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The sample for testing was delivered by the Customer: 09-04-2018
The tests began: 12-04-2018
The tests finished: 20-04-2018

Type of analysis	Method	Results
Microbiological parameters		
Examination of the level of air pollution during the operation of the lamp in a room of 42 m ² and height 2.9 m	Own Methodology Instruction MAS-100 Eco™	[cfu/1 m ³]
- the total number of microorganisms in time 0		547
- the total number of microorganisms after 1 hour		87
- the total number of microorganisms after 2 hours		57
- the total number of microorganisms after 6 hours		17

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-1/1-

Assessment of antibacterial efficacy of UV-C AIR STERILISER ASEPTOR 236

The aim and scope of the research

The aim of the study was to determine the effectiveness of air disinfection by **UV-C AIR STERILISER ASEPTOR 236** (Research report K /166/01/2018) on the basis of the total number of microorganisms examination using aspiration method after 1, 2 and 6 hours steriliser working in a room with an area of 42 m² and height 2.9 m.

Test procedure

The research was conducted in accordance with its own methodology developed at the Laboratory and the manufacturer's instructions MAS-100 ECOTM (Microbiological Air Sampler) in a room with an area of 42 m² and height 2.9 m.

The steriliser was placed in the center of the room and the air pollution was measured 5 meters from the equipment inlet after 1 , 2 and 6 hours of operation. The studies were performed based on the aspiration method using the microbial MAS-100 ECOTM air sampler. Before the lamp was switched on, the general number of microorganisms in the air filling the room was examined. Each time the device was placed on a flat surface, about 90 cm above the floor, facing the outlet up and sucked 1000 liters of air through a perforated plate (time suck approx. 9 minutes). The air stream containing the particles was directed to the surface of PCA agar in a standard Petri dish. After completion of the air sampling cycle, plates were incubated at 30 ° C for 72h. Then grown colonies were counted and the number of microorganisms was determined in 1 m³ air, taking into account the correction of the FELLER statistical conversion table.

Results

Table 1. Level of microbiological air pollution during operation **UV-C AIR STERILISER ASEPTOR 236**

Operating time of the steriliser [h]	The total number of microorganisms [cfu/m ³]
0	547
1	87
2	57
6	17

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